

IN-OFFICE GEL - obj. kód 7343

Aplikácia:

Najprv naneste predbieliaci lak, potom vykonajte 3 aplikácie BEAM In-Office gélu, pričom gél zakaždým nechajte pôsobiť 15 minút.

Tento gél je určený na použitie s lampou BEAM.

Na následnú domácu starostlivosť je určený BEAM WHITENING KIT.

Dôležité informácie na použitie:

Prípravok nie je určený pre osoby mladšie ako 18 rokov.

Výrobok je určený na predaj iba zubným lekárom / stomatologickým odborníkom.

Pri každom cykle použitia musí byť prvá aplikácia vykonaná iba zubným lekárom alebo pod jeho priamym dohľadom, ak je zaistená zodpovedajúca úroveň bezpečnosti. Následne je možné, aby spotrebiteľ výrobok prevzal a dokončil finálne cykly.

Nie je určené pre tehotné alebo dojčiace ženy.

Zabráňte kontaktu s očami. Ak sa výrobok dostane do očí, okamžite ich vypláchnite vodou.

Uchovávať na chladnom a suchom mieste pri teplote 5–25 °C, mimo priameho slnečného žiarenia a mimo dosahu detí.

In-Office Whitening Gel (2,5 ml):

Zloženie: Aqua, Glycerin, Phthalimidoperoxycaproic Acid, Hydrogen Peroxide, Triethanolamine, Carbomer, Sodium Acryloyldimethyltaurate/VP Crosspolymer, Monosodium Citrate, PVM/MA Copolymer, Xanthan Gum, Sodium Saccharin, C12-15 Pareth-5.

Pre-Whitening Varnish (1 ml)

Zloženie: Glycerin, Aqua, Sodium Bicarbonate, Potassium Citrate, Hydroxyapatite, Glycolipids, Sodium Acryloyldimethyltaurate/VP Crosspolymer, Mentha Arvensis Leaf Oil, Stevioside, Mentha Piperita Oil, Potassium Hydroxide, Limonene.